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A newly-developed non-vibrating structure that comprises an energy-saving homogenizing valve.

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Compact sized test equipment for creating samples that makes high-accuracy data collection possible

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### Standard Specifications

- 2-Plungercam Mechanism
- Flow rate 200cc/min
- Pressure 0~50MPa
- Power AC100V 400W
- Pressurization Mechanism  
Automatic pressure control
- Cleaning operation mode

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#### Compact

- Desktop size
- Simple design
- Small number of parts
- Constructed with small, lightweight parts

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#### Highly Operable

- Automatic operations with just 1 button
- Easy breakdown and cleaning
- 100V power allows setup anywhere  
(Please note LABO-02 is AC 200V)
- Output of homo pressure and sample temperature possible

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#### Low Cost

- Tests can be performed with 1 cup
- Compliance and reproducibility with line equipment (large-scale equipment)
- Reduction of work time / laborsaving

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

Desktop homogenizer has enhanced convenience by 100MPa(1,000kg/cm<sup>2</sup>) pressure

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## Standard Specifications

- 2-Plungercam Mechanism
- Flow rate 200cc/min
- Pressure 0~100MPa
- Power AC200V
- Pressurization Mechanism  
Automatic pressure control
- Cleaning operation mode

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### Adoption of a "Cam Method" in the non-vibrating mechanism

Incorporates 2 compact plungers and a non-vibrating mechanism that utilizes cams.

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### Automatic pressure control



After simple operations to set the pressure, the pressure will be maintained automatically.

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
### Cleaning operation mode

The speed of the cleaning water is obtained by adjusting the cleaning water flow rate. Stagnation points in the flow are minimized and care has been taken to make sure that there are no sections that are left unclean.

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
When the fluid volumes pushed out by Plungers 1 & 2 are added, they become constant. A non-vibrating state is achieved in which the pressure and flow rate become constant.

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Supporting the product development of our customers

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Resolves issues that arise when creating samples for research/development.

- 研究開発室にて、お客様のサンプル作製をサポート
- 貸出用ラボ機～生産機をラインナップ
- お客様のニーズに合ったホモジナイザー、ホモバルブをご提案
- At the Research and Development Office, we support our customers' sample creations.
- We have a lineup of loan-out laboratory equipment and production machines for hire
- We offer homogenizers and homogenizing valves that match our customers' needs

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■ 粒度分布計・各種測定器  
Particle size distribution measurement / various measurement devices

■ 各種ホモバルブ  
A variety of homogenizing valves

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**SANMARU MACHINERY CO., LTD.**